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Intended use

The ESR* screw cap tube and commercially available (safety) Luer needles are used together for collecting venous blood. They are used for collecting, transporting and processing venous blood for manual *in vitro* diagnostic determination of the erythrocyte sedimentation rate. The measurement results are obtained according to the Westergren method.

The product is intended for use in a professional environment by qualified medical and laboratory personnel.

*ESR = Erythrocyte sedimentation rate

Product description

The ESR screw cap tube blood collection system consists of a plastic tube with a colour-coded plastic screw cap and a citrate additive. The volume and permissible tolerances of the additive and the ratio of blood to additive comply with the requirements and recommendations of the international DIN EN ISO 6710 standard "Single-use containers for venous blood specimen collection" and the Clinical and Laboratory Standards Institute's Approved Standards (CLS).

Colour code of the ESR screw cap tube blood collection system:

Additive	Letter code	Cap colour based on BS 4851*
ESR screw cap tube	4NC	Violet

*British Standard BS4851:1982 Specification for single use labelled medical specimen containers for haematology and biochemistry, revoked.

SARSTEDT ESR 4NC screw cap tube

The ESR 4NC screw cap tube contains the anticoagulant trisodium citrate and is used for collecting whole blood. Trisodium citrate is present as 0.106 mol/l trisodium citrate equivalent to 3.2% trisodium citrate (3.13% trisodium citrate solution; often rounded up to 3.2%) and amounts to 20% of the nominal volume of an ESR screw cap tube. The mixture ratio of citrate to blood is 1:4 – 1 part by volume citrate and 4 parts by volume blood. Correct filling is absolutely necessary for the analysis. This citrate whole blood is used as specimen material for routine tests, determining the erythrocyte sedimentation rate. Anticoagulation occurs through calcium ion complex formation by citrate.

Safety information and warnings

1. General precautions: Use gloves and other general personal protective equipment to protect yourself from blood and possible exposure to potentially infectious specimen material and transmissible pathogens.
2. Handle all biological specimens and sharp blood collection utensils (needles) according to the guidelines and procedures of your facility. In case of direct contact with biological specimens or of a needlestick injury, consult a doctor due to the risk of transmission of HIV, HCV, or HBV or other infectious diseases. Use (safety) Luer needles that have a built-in needle guard. Always follow the safety guidelines and procedures of your facility.
3. The product is intended for single use. Dispose of all sharp objects (needles) for blood collection in appropriate hazardous biological waste containers.
4. If blood is collected through an intravenous (IV) line, it must be ensured that the line is properly flushed in accordance with facility guidelines (= cleaned of the IV solution) before beginning blood collection. Properly flushing the line prevents anomalous analysis results.
5. Underfilling the ESR screw cap tube blood collection system leads to an incorrect ratio of blood to preparation/additive and can give rise to anomalous analysis results.
6. The blood collected with the ESR screw cap tube blood collection system and then processed is not intended to be re-injected into the human body.
7. Do not use the product after the expiry date. The expiry date of the product is the last day of the month and year indicated.
8. Liquid preservatives and anticoagulants are clear. Do not use if the liquids are cloudy or precipitation is present.

Storage

Store the product at room temperature until use.

Limitations

1. It should be stored at room temperature. The ESR should be determined within the first 4 hours after collecting the blood.
The specimen can be stored in a refrigerator (4 °C) for a longer period (max. 24 hours). The specimen must then be brought to room temperature before being used.
2. Homogenise the blood sample by swirling it carefully (inverting it 5 times) before measuring the erythrocyte sedimentation rate.
3. Measuring requires an ambient temperature of 18–25 °C and protection against vibrations, draughts and direct sunlight must be ensured before measuring.
4. The measuring principle complies with CLSI H2-A5.
5. The measurement accuracy is +/- 1 mm of the measurement path.
6. The mixing ratio of 1:4 has a direct effect on the analysis result and must be adhered to.

Collecting and handling specimens

READ THIS DOCUMENT COMPLETELY BEFORE STARTING VENOUS BLOOD COLLECTION.

Equipment required for specimen collection

1. All required blood collection systems (ESR screw cap tube blood collection system).
2. Commercially available (safety) Luer needle
3. Gloves, gown, goggles or other suitable protective gear to protect against bloodborne pathogens or potentially infectious materials.
4. Labels for patient identification.
5. Disinfectant to clean the collection site (follow the guidelines of the facility for sterile specimen collection for preparing the collection site).
Do not use any alcohol-based cleaning materials if the specimens are to be used for a blood alcohol test.
6. Dry, clean, disposable swabs.
7. A tourniquet.
8. Plaster or dressing.
9. Sharps disposal container for the safe disposal of used material.

Blood collection

1. Apply a venous tourniquet if necessary (max. 1 minute). Prepare the puncture site with a suitable disinfectant. Do not touch the puncture site after disinfection.
Position the patient so that there is easy access to the vein, and the arm or other puncture site on the patient is directed downwards.
2. Open the prepared ESR sample tube, hold it upright and set the cap aside.
3. Puncture the vein with the (safety) Luer needle. Follow the instructions for use for the (safety) Luer needle.
4. Let your specimen material drip into the prepared ESR sample tube. Ensure that the prepared sample tube is not filled with more or less than the nominal volume.
5. Close the prepared ESR sample tube with the corresponding cap. Immediately mix the blood sample thoroughly. Invert the ESR screw cap tube slowly five times.
The air bubble must always move to the other end. Do not shake quickly.
6. Completing blood collection: Withdraw the (safety) Luer needle from the vein and follow the instructions for use for the (safety) Luer needle.
7. Press a clean, dry swab against the puncture site until bleeding stops.
8. As soon as coagulation has set in, you may apply a dressing.
9. Dispose of the (safety) Luer needle in a suitable container for hazardous biological waste.

Performing the ESR measurement

General information

Homogenise the blood sample by swirling it carefully before measuring the erythrocyte sedimentation rate. We recommend using the Sarmix® rotation mixer.

ESR measurement with the ESR screw cap tube

Material:

- Filled ESR screw cap tube
- Sedimentation pipette
- ESR rack
- Online instructions for use: SARSTEDT Manual erythrocyte sedimentation – ESR rack & sedimentation pipettes; GB 562

For determining the ESR, use the SARSTEDT ESR rack with or without a graduated scale on the back. The ESR rack with graduated scale on the back (Art. no. 90.1060) should be used in combination with the sedimentation pipette with a "0" mark (without graduated scale) Art. no.: 86.1996. The ESR rack without back panel (Art. no.: 90.1060.062) should be used in combination with the sedimentation pipette with graduated scale (Art. no.: 86.1996.062). When determining the ESR, follow the instructions for use of the respective SARSTEDT ESR rack & sedimentation pipettes (GB 562; SARSTEDT Manual erythrocyte sedimentation – ESR rack & sedimentation pipettes at: www.sarstedt.com/en/download/instructions-for-use).

Disposal

1. The general hygiene guidelines and regulations for the proper disposal of infectious material must be observed and complied with.
2. Disposable gloves prevent the risk of infection.
3. Contaminated or filled blood collection systems must be disposed of in suitable containers for hazardous biological waste that can be subsequently autoclaved and incinerated.
4. Dispose of potentially contaminated material according to the regulations and guidelines of the facility.

Product-specific standards and guidelines as amended / updated

CLSI* GP39 "Tubes and Additives for Venous and Capillary Blood Specimen Collection", Approved Standard.

CLSI* GP41 "Collection of Diagnostic Venous Blood Specimens", Approved Standard.

CLSI* GP44 "Procedures for the Handling and Processing of Blood Specimens for Common Laboratory Tests", Approved Guideline.

CLSI* H2-A5 2011: Procedures for the Erythrocyte Sedimentation Rate Test, 5th Edition Item Order Code: H02ASE.

WHO/DIL/LAB/99.1 Rev02 "WORLD HEALTH ORGANIZATION, et al. Use of anticoagulants in diagnostic laboratory investigations. Geneva: World Health Organization, 2002".

Standard for containers for venous blood collection:

DIN EN ISO 6710: Single-use containers for human venous blood specimen collection

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ICSH recommendations for measurement of erythrocyte sedimentation rate. J. Clin. Pathol. 1993 46:198-203

*CLSI (Clinical and Laboratory Standards Institute)

Key for symbols and labels:

Article number

Batch number

Use by

CE mark

In vitro diagnostic device

Follow the instructions for use

If reused: Risk of contamination

Keep away from sunlight

Store in a dry place

Manufacturer

Country of manufacture

Technical modifications reserved.

All serious incidents relating to the product shall be notified to the manufacturer and the competent national authority.

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*British Standard BS4851:1982 Specification for single use labelled medical specimen containers for haematology and biochemistry, изтетен.

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CLSI* GP39 "Tubes and Additives for Venous and Capillary Blood Specimen Collection", Approved Standard.
 CLSI* GP41 "Collection of Diagnostic Venous Blood Specimens", Approved Standard.
 CLSI* GP44 "Procedures for the Handling and Processing of Blood Specimens for Common Laboratory Tests", Approved Guideline.
 CLSI* H2-A5 2011: Procedures for the Erythrocyte Sedimentation Rate Test, 5th Edition Item Order Code: H02ASE.
 WHO/DIL/LAB/99.1 Rev02 "WORLD HEALTH ORGANIZATION, et al. Use of anticoagulants in diagnostic laboratory investigations. Geneva: World Health Organization, 2002".

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CLSI* GP39 " Tubes and Additives for Venous and Capillary Blood Specimen Collection" , Approved Standard.

CLSI* GP41 " Collection of Diagnostic Venous Blood Specimens" , Approved Standard.

CLSI* GP44 " Procedures for the Handling and Processing of Blood Specimens for Common Laboratory Tests" , Approved Guideline.

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WHO/DIL/LAB/99.1 Rev02 " WORLD HEALTH ORGANIZATION, et al. Use of anticoagulants in diagnostic laboratory investigations. Geneva: World Health Organization, 2002" .

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