



Protein Residue Test

General Protein

AllerSnap®

High-sensitivity broad spectrum protein

🕒 15 min

PRO-Clean

Protein

🕒 1 - 10 min

SpotCheck Plus

Glucose & Lactose

🕒 1 min

Lateral Flow Tests

AlerTox® Sticks

Trace allergen contamination in food, liquid, and surfaces

Allergens that can be detected with AlerTox Sticks:

- 🥜 Almond
- 🥛 Casein
- 🦀 Crustacean
- 🥚 Egg
- 🐟 Fish
- 🌰 Hazelnut
- 🥛 β - Lactoglobulin
- 🌱 Mustard Seed
- 🥜 Peanut
- 🌱 Sesame Seed
- 🥛 Soy
- 🥛 Total Milk
- 🌰 Walnut

🕒 10 min

GlutenTox® Pro

AOAC-certified gluten detection harnessing the power of the G12 antibody



🕒 10 min

GlutenTox® Sticks Plus

Detection limit 3 ppm (< 20 ppm is considered "gluten-free")

🕒 10 min

Hygiena® Cube

Readable and trackable quantitative or Pass/Fail results

Allerflow | Gluten

High-sensitivity gluten specific

🕒 10 min

ELISA

AlerTox® ELISA

🥜 AlerTox ELISA Almond	🥚 AlerTox ELISA Lysozyme
🥛 AlerTox ELISA BLG	🌰 AlerTox ELISA Macadamia
🥛 AlerTox ELISA Casein	🥛 AlerTox ELISA Milk
🌱 AlerTox ELISA Cashew	🌱 AlerTox ELISA Mustard
🥥 AlerTox ELISA Coconut	🥚 AlerTox ELISA Ovalbumin
🦀 AlerTox ELISA Crustacean	🥜 AlerTox ELISA Peanut
🥚 AlerTox ELISA Egg	🌱 AlerTox ELISA Pistachio
🐟 AlerTox ELISA Fish	🌱 AlerTox ELISA Sesame
🌰 AlerTox ELISA Hazelnut	🥛 AlerTox ELISA Soy
🐾 AlerTox ELISA Lupine	🌰 AlerTox ELISA Walnut

🕒 60 min

GlutenTox® ELISA

Gluten Tox® ELISA Rapid G12
GlutenTox® ELISA Competitive G12

🕒 90 min

PCR

foodproof®

🌱 Celery

🌾 Gluten

🌰 Hazelnut

🥜 Peanut

🥛 Soya

🌰 Walnut

🕒 180 min