



ELECTROSTATIC PITH BALLS

Cat: EM1775-001 (with conductive coating) EM1774-001 (plain)

DESCRIPTION:

A pack of 12x electrostatic 'Pith Balls' for performing electrostatic experiments.

These pith pieces are approximately 'ball' shaped and are available either as plain pith or treated with a conductive coating.

EM1775-001 conductive

EM1774-001 non conductive



Physical size: approx 10mm diameter

Pack/12 in plastic bag

Used for general electrostatic experiments. The conductive 'balls' carry the electrostatic charge all over their surface and the whole surface charges or discharges.

Manufactured in Australia