



## PROOF PLANE – flat disc or sphere

EM2520-001 20mmD disc on handle EM2521-001 (sphere on handle)

### DESCRIPTION:

The classic Proof Plane is a 20mm diameter metal disc or metal sphere attached to a insulating handle. It is used to collect electrostatic charges from devices and to deposit them on other devices.

Design can be either brass or aluminium alloy and the handle can be either acrylic or polystyrene. Performance is exactly the same.

**EM2521-001 & EM2520-001 proof plane / sphere**



Physical size: handle: 10x150mm DxL 20mmD disc or sphere Weight: 0.05 kg

Designed and manufactured in Australia