


**COLOUR FILTER SET - 10 colours, 100x100mm****Cat: HL1882-001 set/10 colours, 100x100mm plain, unmounted****DESCRIPTION:**

The IEC set of 10x colour filters is used in general optics and in the study of colour mixing. Material is cellulose acetate, as used for lighting in stage productions.

**HL1882-001 set / 10 colour filters****Physical size: 100x100mm LxW****Weight: 0.02 kg / set**

Code	Colour	Size
HL1882-002	Primary red	100x100mm
HL1882-003	Primary green	100x100mm
HL1882-004	Magenta	100x100mm
HL1882-005	Violet	100x100mm
HL1882-006	Primary blue	100x100mm
HL1882-007	Yellow	100x100mm
HL1882-008	Cyan	100x100mm
HL1882-009	Orange	100x100mm
HL1882-010	Bright blue	100x100mm
HL1882-011	Yellow chrome	100x100mm

**Designed and manufactured in Australia**