



## POLAROID FILTERS - pair

**Cat: HL1883-010 (50mm sq) HL1883-020 (75mm sq) HL1883-030 (100mm sq)**


### DESCRIPTION:

The IEC Linear Polarising films are available in pairs of squares of various sizes. They can be used as squares or they can be cut into shapes or strips to suit various experiments.

The transmission is 80% when polaroids are aligned. When misaligned at 90°, the transmission reduces to 0.2%.

**NOTE:** HL1883-010, HL1883-030, and HL1883-030 are pairs of filters cut to size. HL1883-040 is a single sheet of polaroid material that can be cut to the shapes desired.

**HL1883-010    HL1883-020    HL1883-030    HL1883-040**



Code	Type		Size & thickness
HL1883-010	Linear polariser	pair	50x50mm x 0.18mm thick
HL1883-020	Linear polariser	pair	75x75mm x 0.18mm thick
HL1883-030	Linear polariser	pair	100x100mm x 0.18mm thick
HL1883-040	Linear polariser	Single sheet or pieces	Totalising 300x300mm x 0.18mm thick

**Industrial Equipment & Control Pty. Ltd.**