

## **LIGHT SOURCE – low power**

**HL2011-001 (MES on base w/4mm sockets)**

### **DESCRIPTION:**

IEC manufactures several low power light sources for general basic electricity experiments. These are low cost and simple lamp holders with the lamp exposed.

This item is the spare part from the IEC Electricity Kit (EM1763-001).

See also HL2010-001, HL2012-001 and HL2012-101 for the types of low power light source holders where the lamps are exposed or are inside a simple vial housing.

**HL2011-001** holds an MES lamp on a large stable base using a printed circuit board to permanently connect the 4mm sockets to the lamp socket.

**HL2011-001 (high quality MES holder)**



**Physical size: 75x75x25mm LxWxH**

**Weight: 0.06 kg**

**Designed and manufactured in Australia**