




SCREEN – frosted acrylic

Cat: HL3740-001 200x200mm, with foot.

DESCRIPTION:

This small acrylic screen is for general optics in the classroom. The ‘frosting’ on one side is a fine sandblast that provides a good image. The simple foot is formed by a strong pin fixed through the acrylic to rest on the table on each side of the screen.

HL3740-001 acrylic screen, frosted



Physical size: 200x200x3mm LxWxTh

Weight: 0.15 kg

Fine frosted surface on one side used for optics and lenses and the the projection of images on a screen.

Designed and manufactured in Australia