



PENDULUM BOBS - set/3

Cat: MF2290-001 set of 3 brass pendulum bobs

DESCRIPTION:

The IEC **Set of Brass Pendulum Bobs** consists of 3x Bobs:

1x Brass bob: 25mm diameter

1x Brass bob: 20mm diameter

1x Brass bob: 12mm diameter

These diameters are nominal only and, if experiments involving mass are performed, the actual diameter of each bob should be measured accurately. Each bob has a 1.5mm diameter hole drilled right through the diameter with one side made larger so that a knot in a cord can be pulled into the larger hole to hold firmly and also be out of sight.

NOTE: The picture below shows also a steel and a plastic bob, but these are not included in the set/3 brass bobs.

MF2290-001 (set/3 brass bobs)



Physical size: Small plastic bag containing 3x Bobs **Weight:** 0.21 kg

Designed and manufactured in Australia