



PULLEY IN SHEAVE - single, double, triple or quad.

Cat: MF2560-001 MF2570-001 MF2580-001 MF2582-001

DESCRIPTION:

The IEC 'Pulley in Sheave' is a high quality and strong nylon pulley free running in a metal sheave. The pivot shaft is silver steel and the 50mm diameter pulley is moulded from low mass nylon. Friction is very low so that experiments remain accurate.

A very strong hook at each end permits the use of heavy weights without causing damage or distortion. The pulleys can be dismantled if any repair or changes to sheaves is required.

MF2560-001 single / MF2570 double / MF2580 triple / MF2582 quad



Physical size: each pulley: 50x12mm DxTh

Weight: 0.03 kg/sheave

IEC manufactures a range of sheaved pulleys from single pulley to four pulley models.

- MF2560-001 Single sheave
- MF2570-001 Double sheave
- MF2580-001 Triple sheave
- MF2582-001 Quad sheave

Other IEC pulleys include bench mounting, rod mounting.

Designed and manufactured in Australia