



## **PULLEY – bench mounting**

**Cat: MF2590-001 bench mounting pulley**

### **DESCRIPTION:**

The IEC 'Pulley on Clamp' is a high quality and strong nylon pulley free running on a small metal 'G' clamp. The clamp can open to a maximum of 25mm. The 50mm diameter pulley runs parallel to the bench top.

**NOTE:** For larger bench clamps refer to LB0858-001. This item can open to fit to very thick benches and can accept any type of rod to be clamped, including retort rods or pulley on rod MF2600-020 and MF2600-001.

**MF2590-001 small bench pulley**



**Physical size: 50x75mm DxL**

**Weight: 0.04 kg**

**Designed and manufactured in Australia**