

**WEIGHTS – slotted, brass, 50g, 25g, 20g**

MF4300-002 50g, MF4300-005 25g, MF4300-006 20g,

individual slotted brass weights**DESCRIPTION:**

The IEC ‘**Slotted Brass Weights**’ are solid brass discs with a shaped slot that permits captive engagement on a 5mm diameter rod in the form of a ‘carrier’ with a strong hook fitted permanently into the top. Each weight is marked 50g, 25g or 20g.

Each weight and the carrier is calibrated to an accuracy is better than 1%.

Several types of slotted weights are manufactured by IEC.

MF4300-002 50g, MF4300-005 25g, MF4300-006 20g**Physical size: 30mmDx various thickness Weights: 0.05 kg, 0.025 kg, 0.02 kg****Designed and manufactured in Australia**