



WEIGHTS – slotted, brass, set/9 x 100g on carrier

MF4301-001 slotted brass weights, set/9 100g on carrier

DESCRIPTION:

The IEC ‘**Slotted Brass Weights**’ are solid brass discs with a shaped slot that permits captive engagement on a 5mm diameter rod in the form of a ‘carrier’ with a strong hook fitted permanently into the top. The set of 9 weights fitted to a 100g carrier makes a total weight of 1000g. Each weight in the set is marked 100g together with the brass carrier with hook.

Each weight and the carrier is calibrated to an accuracy is better than 1%.

Several types of slotted weights are manufactured by IEC.

MF4301-001 set / 9x 100g on carrier



Physical size: 30x217mm DxH

Weight: 1.0 kg

Designed and manufactured in Australia