



WEIGHTS – slotted, brass, set/5

MF4304-001 slotted brass weights, set/5: 10g, 2x20g, 50g, 100g

DESCRIPTION:

The IEC ‘**Slotted Brass Weights**’ are solid brass discs with a shaped slot that permits captive engagement on a 5mm diameter rod in the form of a ‘carrier’ with a strong hook fitted permanently into the top. The set of 5 weights does not include a carrier.

The set is: 1x 10g, 2x 20g, 1x 50g, 1x 100g to make total weight of 200g.

Each weight is calibrated to an accuracy is better than 1%.

Several types of slotted weights are manufactured by IEC.

MF4304-001 set / 5x slotted weights



Physical size: 30mmDx various thickness Weight: 0.20 kg

Designed and manufactured in Australia