



## **WEIGHTS – slotted, brass, set/3**

**MF4304-101 slotted brass weights, set/3: 5g, 10g, 1x20g**

### **DESCRIPTION:**

The IEC ‘**Slotted Brass Weights**’ are solid brass discs with a shaped slot that permits captive engagement on a 5mm diameter rod in the form of a ‘carrier’ with a strong hook fitted permanently into the top. The set of 3 weights does not include a carrier.

The set is: 1x 5g, 1x 10g, 1x 20g to make total weight of 35g.

Each weight is calibrated to an accuracy is better than 1%.

Several types of slotted weights are manufactured by IEC.

### **MF4304-101 set / 3x slotted weights**



**Physical size: 30mmDx various thickness**

**Weight: 0.03 kg**

**Designed and manufactured in Australia**