

# CULTURA<sup>®</sup>M



CULTURA<sup>®</sup>M

## The CULTURA<sup>®</sup>M

Small incubator is suitable for the incubation of staining or immersion objective slides such as Uricult, Hygicult, Easicult, (all products of Orion Diagnostica, Finland) or regular Petri dishes for the evidence of bacteria, fungi and yeast.

The CULTURA<sup>®</sup>M has been used for more than 20 years in the fields of medical and dental diagnostics, food bacteriology and hygiene monitoring, monitoring of coolants lubricants and industrial liquids as well as of drinking water, industrial process and swimming pool water.

The CULTURA<sup>®</sup>M is user-friendly, tried-and-tested, durable and reliable.

## Technical Data

- Trademark: CULTURA<sup>®</sup>
- Part numbers: 70700 – 70709, 70711, 70713, 70714, 70700R – 70709R, 70711R, 70713R, 70714R
- Temperature range: +25 °C bis +45 °C
- Temperature stability: +/- 1°C
- Voltage/power output  
220-240 V / 30 W / 0.2 A  
100-220 V / 30 W / 0.2 A  
12 V / 20 W / 0.2 A
- Frequency: 50/60 Hz
- Protection class: II
- Case: hardly inflammable PC & ABS
- Door: transparent PC
- Thermometer: up to +60 °C (in special tube)
- Base cover or Multirack: ABS
- Warranty: two year from initial installation date
- External dimensions: L x B x H: 310 x 168 x 155 mm
- Internal dimensions: L x B x H: 220 x 150 mm x 120
- Weight: 1,1 kg