

Tube Roller Mixer

Features & Benefits

- Tube roller
- Rolling and rocking mixing action for complete mix
- Can be used in cold rooms and incubators
- DC brushless motor



400.104.005

Tube Roller Mixer (tubes not supplied).

Technical Specifications:

Product Code	400.104.005
Voltage	240V/60Hz
Power	25W
Motor Type	DC Brushless
Mixing Action	Rolling and Rocking
Amplitude	24mm
Maximum Load Capacity	4kg
Number of Rollers	6
Roller Length	280mm
Speed Range	0 to 70rpm
Speed Display	Scale
Operation	Continuous
Dimensions	450 x 260 x 120mm
Weight	4.5kg
Permissible Ambient Temperature	5 to 40°C
Permissible Relative Humidity	80%
Protection Class	IP21

Ordering Information:

Product Code	Description	Pack Size
400.104.005	Tube Roller	Each



400.104.005 Control Panel



400.104.005 Loaded with Labco 650.550.015 (4 x 15mL) & 650.550.050 (6 x 50mL) Centrifuge Tubes (tubes not supplied).