

Autoclave Tape



Autoclave tape is a special type of adhesive tape primarily used in laboratory or medical settings to indicate whether an item or package has successfully gone through the sterilisation process in an autoclave.

FEATURES + BENEFITS

Visual Sterilisation Verification

The heat sensitive chemical indicators change colour when exposed to high temperature and pressure conditions inside an autoclave. Black stripes will appear on the tape to indicate the sterilisation cycle was successful.

Preventing Cross Contamination

Can prevent the risk of cross contamination in healthcare settings, by ensuring medical instruments are sterile before use.


Cost Effective

Compared to other sterilisation indicators, autoclave tape is a cost effective solution that's simple and easy to use.

Lead Free



Three sizes to choose from



24mm

18mm

12mm

Ordering Information

CODE	DESCRIPTION	PACK SIZE
500.750.012	Tape Autoclave 12mm x 50m	Roll
500.750.018	Tape Autoclave 18mm x 55m	Roll
500.750.024	Tape Autoclave 24mm x 55m	Roll

