

Plasticware

Pasteur (Transfer) Pipettes

Features & Benefits

- Manufactured from non-toxic, Low-Density Polyethylene (LDPE)
- Ideal for transferring and dispensing liquids
- Eliminates the risk of cross contamination
- Excellent transparency
- Uniform wall thickness and precise graduations
- Choice of packaging; bulk, sterile in inner packs of 20, sterile individually wrapped



Technical Specifications:

Capacity	5mL
Total Length	160mm
Graduated	Up to 1mL
Graduation Increments	0.25mL
Drops per 1mL	25
Stem Diameter	5mm
Bulb Draw	3.5mL



650.050.501
Bulk
Non-sterile



650.050.511
Inner Packs of 20
Sterile




650.050.521
Individually Wrapped
Sterile

Ordering Information:

Product Code	Packaging	Sterile	Pack Size
650.050.501	Bulk	No	PK500
650.050.511	Inner packs of 20	Yes	PK500
650.050.521	Individually wrapped	Yes	PK500

Features & Benefits

- Manufactured from non-toxic, Low-Density Polyethylene (LDPE)
- Ideal for transferring and dispensing liquids
- Eliminates the risk of cross contamination
- Excellent transparency
- Uniform wall thickness and precise graduations
- Choice of packaging; bulk, sterile in inner packs of 20, sterile individually wrapped



Technical Specifications:

Capacity	7mL
Total Length	160mm
Graduated	Up to 3mL
Graduation Increments	0.5mL
Drops per 1mL	21
Stem Diameter	7.8mm
Bulb Draw	3.5mL



650.050.503
Bulk
Non-sterile



650.050.513
Inner Packs of 20
Sterile



650.050.523
Individually Wrapped
Sterile



Ordering Information:

Product Code	Packaging	Sterile	Pack Size
650.050.503	Bulk	No	PK500
650.050.513	Inner packs of 20	Yes	PK500
650.050.523	Individually wrapped	Yes	PK500