



LABORATORY OVENS



BENCHTOP LABORATORY OVENS



(AMBIENT +10°C TO +200°C)

(Image: TO-30G)



DIGITAL LABORATORY OVENS

Lab Ovens with Digital Temperature Control, Max +200°C

- Stainless steel internal liner (easy clean surface and corrosive resistant)
- Fibreglass insulation for energy efficiency
- Gravity Convection (No Fan)
- Stainless steel sheathed electric heating element
- Digital set point controller featuring error codes in the event of malfunction
- Monitoring port hole
- Removable and Adjustable chrome plated shelving
- Australian Made



PRODUCT DETAILS

These laboratory ovens are designed and manufactured in Australia by Thermoline for all laboratory and light industrial functions. Both models feature a stainless steel internal liner which is easy to clean and corrosive resistant. A protected calrod heating element is located below a false floor providing a gentle, natural convection heat throughout the workspace. Fibreglass insulation is in all walls for power efficiency and a digital PID controller featuring error codes in the event of malfunction. All Thermoline ovens are very well engineered making them reliable and long lasting.

For larger models, please view Thermoline's 'Premium' Dehydrator range.

Optional Accessories

- Additional shelving available
- + Additional monitoring cable port holes
- Controlled temperature ramp to setpoint

Safety Features

- ! Auto overtemperature safety (Failsafe)
- Over current protection

SPECIFICATIONS

| Model | TO-30G | TO-150G |
|--------------------------------|-------------------------|-------------|
| Internal Dimensions (WxDxH mm) | 300x285x420 | 500x510x600 |
| External Dimensions (WxDxH mm) | 420x400x735 | 620x670x930 |
| Capacity (litres) | 30 | 150 |
| Shelves | 3 | 3 |
| Temperature Range | Ambient +10°C to +200°C | |
| Electrical (W) | 500 | 2000 |
| Weight (kg) | 35 | 78 |

Thermoline Scientific have been manufacturing and distributing high quality laboratory and scientific testing equipment since 1970. Over this time, Thermoline has grown to be a leading brand in the science industry, with our equipment being used in small and large Hospitals, Universities and Research Laboratories across Australia and the Asia Pacific region.

We're proud to say that we are 100% Australian owned and operated.

Head Office Phone (02) 9604 3911

Head Office Fax (02) 9725 1706

Email hello@thermoline.com.au

Web www.thermoline.com.au





Thermoline Scientific have been manufacturing and distributing high quality laboratory and scientific testing equipment since 1970. Over this time, Thermoline has grown to be a leading brand in the science industry, with our equipment being used in small and large Hospitals, Universities and Research Laboratories across Australia and the Asia Pacific region.

We're proud to say that we are 100% Australian owned and operated.