

Product Specification

Product Name: NUTRIENT AGAR LR

Alternate Name(s) Non-animal origin

Description

Light yellow to yellow in appearance, may have a slight green tinge

Properties

Chemical Formula:

Molecular Weight:

General Information:

Product Code: NL023

CAS No.

Quality Specification

Assay:

Specific Properties and Impurities [Typical levels]:

Gelling	Firm, compared with 1.5% agar gel
pH @ 25°C	7.4 ± 0.2
Cultural response at 35-37 deg. C for 18-48 hr:	
Organism:	
E. coli ATCC 25922	Luxuriant
P. aeruginosa ATCC 27853	Luxuriant
S. aureus ATCC 25923	Luxuriant
S. typhi ATCC 6539	Luxuriant
S. pyrogenes ATCC 19615	Luxuriant

Hazard and Safety Data

UN Group:	None Allocated
Class:	None Allocated
UN Number:	None Allocated
Hazchem code:	None Allocated
CS MSDS Code:	
Poison schedule:	Not scheduled
Emergency Procedure Guide No.:	None Allocated

Chem-Supply Pty Ltd - An ISO 9001 Accredited Company

**38 - 50 Bedford Street, Gillman SA 5013, Australia ABN 19 008 264 211 PO Box 201, Port Adelaide SA 5015, Australia
Telephone +61 8 8440 2000 Fax +61 8 8440 2001 E-mail: sales@chemsupply.com.au Web: www.chemsupply.com.au**

Chem-Supply does not warrant that this product is suitable for any use or purpose. The user must ascertain the suitability of the product for any intended purpose. Preliminary testing of the product before use or application is recommended. Any reliance or purported reliance upon Chem-Supply Pty Ltd with respect to any skill or judgement or advice in relation to the suitability of this product for any purpose is disclaimed. Except to the extent prohibited at law, any condition implied by any statute as to the merchantable quality of this product or fitness for any purpose is hereby excluded. This product is not sold by description. Where the provisions of Part V, Division 2 of the Trade Practices Act apply, the liability of Chem-Supply Pty Ltd is limited to the replacement or supply of equivalent goods or payment of the cost of replacing the goods or acquiring equivalent goods.